



## MAGPRO R30

Biphasic, Biphasic Burst (Theta Burst)

An advanced high-performance TMS magnetic stimulator designed primarily for clinical use as well as research. The MagPro R30 performs repetitive transcranial magnetic stimulation and runs many complex protocols successively; up to 20,000 pulses in one session.

With the add-on Theta Burst option, it is possible to run theta burst stimulation. An add-on 60 pps option is also available.

### Contact

Website: [www.magventure.com](http://www.magventure.com)  
Email: [info@magventure.com](mailto:info@magventure.com)

**Max repetition rate**  
30 pps, 60 pps, 80 pps

**Pulse mode**  
Standard

**Waveform**  
Biphasic, Biphasic Burst (Theta Burst)

**Current direction**  
Normal

**Ramp up**  
Yes

### Intended use

**Diagnostic Use** The MagPro R30 is intended as an electro-physiological aid to assess diagnosis and to monitor diseases of the central and peripheral nervous system, based on the use of Motor Evoked Potentials (MEP). **Therapeutic Use** Treatment of Major Depressive Disorder in adult patients who have failed to achieve satisfactory improvement from two prior antidepressant medications, at or above the minimal effective dose and duration in the current episode.

### Mechanical Data

HxWxD	210x530x400 mm/ 8.3x20.9x15.7 in
Weight	36 kg/79 lbs

### Environmental Data

Operating temperature	10-30°C/50-85°F
Storage temperature	5-50°C/41-122°F

Operating humidity	30-60% RH
Storage humidity	20-80% RH

## Power Supply

Mains voltage	230V~,50/60 Hz according to IEC 60601-1
Maximum power consumption	2700 VA, operation from 100 V~ or 120 V~ through transformer

Ordering number: [9016E0721](#)